



How much DNA do I need to send for my genotyping project?

DNA can be shipped to LGC in a dried-down or liquid format. Please refer to the relevant section below depending on how you would prefer to ship your samples:

Section 1: Dry DNA samples

Section 2: Wet DNA samples

The values stated here are intended as a guideline, and will vary based on the quality of your DNA samples.

When you send in your DNA samples, please include as much information as possible i.e. DNA concentration, method used for extraction, method used for DNA quantification.

Please note: it is always preferable to send in DNA at a higher mass or concentration than the minimum requirement. All DNA will be tested in-house before commencing your genotyping project to determine the most suitable dilution. Any unused DNA can be stored at LGC in the short term for future genotyping projects, or arrangements can be made to return surplus DNA to the customer.

Please ensure that the concentration of your DNA samples is uniform. LGC are able to normalise your samples but this would incur an additional charge.

A range of methods can be used for quantification of DNA samples, including spectrophotometry and PicoGreen. Spectrophotometers have a tendency to overestimate the quantity of DNA present whilst PicoGreen is more exact. All DNA quantities / concentrations in this document assume that samples have been quantified using spectrophotometry.

Please note: for DNA that has been extracted from saliva / buccal swabs, bacterial DNA will also be present and will affect the accuracy of quantification data.

1. Dry DNA samples

The mass of DNA required per sample per SNP will differ depending on the genome size of your study organism. If you are working with human DNA (genome size ~3000 Mbp) or a species with a genome size in the range of 2000 – 3500 Mbp, we require 10 ng of good quality DNA per sample per SNP (except for DNA quantified using PicoGreen, where a minimum of 5 ng per sample per SNP is acceptable).

Please see the table below for a general guide to the amount of DNA required when sending in dried down samples. The exact amounts required will vary based on the starting material and the extraction procedures used.

When calculating the mass of DNA required for dried down samples, based on the number of SNPs to be run, always send an additional 20 - 25% mass. This will account for potential losses during re-suspension of samples.

Genome size	Minimum mass of DNA required per sample per SNP	Example species
100 Mbp – 750 Mbp	2 ng (N.B. this is the minimum mass that we accept)	<i>Oryza sativa</i> (rice), <i>Arabidopsis thaliana</i> , <i>Drosophila melanogaster</i> (fruit fly)
750 Mbp – 2000 Mbp	2.5 – 6.5 ng	<i>Solanum tuberosum</i> (potato), <i>Glycine max</i> (soybean)
2000 Mbp – 3500 Mbp	10 ng	<i>Homo sapiens</i> (human), <i>Zea mays</i> (maize/corn), <i>Helianthus annuus</i> (sunflower), <i>Mus musculus</i> (house mouse), <i>Danio rerio</i> (zebrafish)
3500 Mbp – 5000 Mbp	11 – 16 ng	<i>Lens culinaris</i> (lentil)
5000 Mbp – 10000 Mbp	16 – 33 ng	
10000 Mbp – 15000 Mbp	33 – 50 ng	
15000 Mbp – 20000 Mbp	50 – 66 ng	<i>Allium cepa</i> (onion)
20000 Mbp – 30000 Mbp	66 – 100 ng	<i>Pinus cembra</i> (pine)
30000 Mbp – 50000 Mbp	100 – 166 ng	<i>Hordeum vulgare</i> (barley)



2. Wet DNA samples

The concentration of DNA required will vary based upon the genome size of your study organism. The minimum DNA concentration that we accept is 5 ng / μ L (except for DNA quantified using PicoGreen, where a minimum of 3.3 ng / μ L is acceptable).

Please note: it is always advisable to send DNA at a higher concentration than the minimum requirement. This allows for a range of dilutions to be tested in-house and reduces the requirement for dead volumes. Any unused DNA can be stored at LGC in the short term for future genotyping future genotyping projects, or arrangements can be made to return surplus DNA to the customer.

Genome size	Minimum concentration of DNA required	Example species
100 Mbp – 750 Mbp	5 ng / μ L	<i>Oryza sativa</i> (rice), <i>Arabidopsis thaliana</i> , <i>Drosophila melanogaster</i> (fruit fly)
750 Mbp – 2000 Mbp	5 ng / μ L	<i>Solanum tuberosum</i> (potato), <i>Glycine max</i> (soybean)
2000 Mbp – 3500 Mbp	5 ng / μ L	<i>Homo sapiens</i> (human), <i>Zea mays</i> (maize/corn), <i>Helianthus annuus</i> (sunflower), <i>Mus musculus</i> (house mouse), <i>Danio rerio</i> (zebrafish)
3500 Mbp – 5000 Mbp	5 – 11 ng / μ L	<i>Lens culinaris</i> (lentil)
5000 Mbp – 10000 Mbp	11 – 22 ng / μ L	
10000 Mbp – 15000 Mbp	22 – 33 ng / μ L	
15000 Mbp – 20000 Mbp	33 – 44 ng / μ L	<i>Allium cepa</i> (onion)
20000 Mbp – 30000 Mbp	44 – 67 ng / μ L	<i>Pinus cembra</i> (pine)
30000 Mbp – 50000 Mbp	67 – 111 ng / μ L	<i>Hordeum vulgare</i> (barley)

The table below gives a general guide to the concentration of DNA required, based on genome size of the study organism. LGC requires 1.5 μ L DNA per sample per SNP at an appropriate concentration for the genome size of your study organism. This volume is calculated based on the minimum concentration; if higher concentrations are supplied, the volume requirement is reduced.

Due to laboratory handling procedures, we do not accept volumes of less than 20 μ L. If you only have a small number of SNPs to run, you will still need to supply 20 μ L of starting material. We also require a dead volume; this is an extra volume of DNA required for the robotic pipetting process that we use. This volume will vary based on the total number of SNPs in your genotyping project. Please see the table below for an outline of the volumes of DNA to send.

To calculate the volume of DNA required (assuming DNA is at the minimum required concentration):

Number of SNPs x 1.5 μ L + dead volume = Total volume of DNA required

Number of SNP SNPs to be run on your samples	Volume of wet DNA (μ L) to be sent (exclusive of 'dead volumes')	Amount of additional DNA 'dead volume' required (μ L) for robotics	Total volume of DNA required (μ L)*
Up to 5	10	10	20
Up to 10	15	10	25
Up to 25	38	20	58
Up to 50	75	20	95
Up to 100	150	30	180
Up to 200	300	30	330

*N.B. Ensure that the DNA is at or above the minimum concentration required for your study organism.

Please note: If a large volume of DNA is required for your desired number of SNPs, please send this DNA at a higher concentration than the minimum acceptable concentration to reduce the volume required. This will enable the whole sample to be sent in one plate (maximum well volume 0.8 mL). If you are supplying your DNA in greater concentrations than required for your number of SNPs, the DNA samples can be diluted in house. In these instances, the in-house dilution reduces the need for dead volumes.

Calculating the mass of DNA required according to genome size

If your study organism has a:

- **Larger genome size** – you will need to send a greater quantity of DNA. To calculate this, divide the genome size of your organism by the size of the human genome (3000 Mbp), and use the resulting number to multiply the amount of DNA that you need to send.
e.g. *Triticum aestivum* (wheat): 15966 Mbp
 $15966 \text{ Mbp} / 3000 \text{ Mbp} = 5.3$
You will need 5.3 times as much wheat DNA per sample per SNP = $10 \text{ ng dry DNA} \times 5.3 = 53 \text{ ng}$
- **Smaller genome size** – you can send a lower quantity of DNA. To calculate this, divide the human genome size (3000 Mbp) by the genome size of your organism and use the resulting number to divide the amount of DNA that you need to send.
e.g. *Oryza sativa* (rice): 441 Mbp
 $3000 \text{ Mbp} / 441 \text{ Mbp} = 6.8$
You will need 6.8 times less rice DNA per sample per SNP = $10 \text{ ng dry DNA} / 6.8 = 1.5 \text{ ng}$

Please note, we do not accept below 2 ng DNA for any species.

Whole Genome Amplification

If you do not have sufficient genomic DNA for the number of SNPs that you wish to run, it is possible to perform Whole Genome Amplification (WGA). WGA is a PCR technique that is used to produce large quantities of DNA from a small amount of starting material. There are a number of methods for WGA, and LGC favours use of the primer extension pre-amplification (PEP) technique. PEP employs the use of randomly synthesised 15-mer oligonucleotides, referred to as polyN15, that bind at sites throughout the genome and act as primers to enable DNA replication. LGC uses our in-house enzyme, Klear Taq, and buffer system to perform WGA reactions.

To perform WGA on your samples, LGC requires a minimum of 50 ng of genomic DNA per sample. If the quality of the starting material is good, the product is typically amplified to a concentration of 500 - 1000 times higher than that of the starting material.

Quantification of DNA samples using PicoGreen

If you have used PicoGreen for quantification of your DNA samples, you may be able to send a lower concentration of DNA to LGC for your genotyping project. This is because PicoGreen quantification is more exact than spectrophotometry.

* For dry DNA samples quantified by PicoGreen, work on the basis of 5 ng per sample per reaction rather than 10 ng per sample per reaction

*For wet DNA samples quantified by PicoGreen, work on the basis of 3.3 ng / μL per reaction rather than 5 ng / μL per sample per reaction.

www.lgcgroup.com/genomics
genomics@lgcgroup.com

Germany

LGC
Ostendstr. 25 • TGS Haus 8
12459 Berlin • Germany

Tel: +49 (0)30 5304 2200
Fax: +49 (0)30 5304 2201

United Kingdom

LGC
Unit 1-2 Trident Industrial Estate • Pindar Road
Hoddesdon • Herts • EN11 0WZ • UK

Tel: +44 (0)1992 470 757
Fax: +44 (0)1438 900 670

United States

LGC
100 Cummings Center • Suite 420H
Beverly • MA 01915 • USA

Tel: +1 (978) 232 9430
Fax: +1 (978) 232 9435