



# RHH Liquid Nitrogen Store: Operation & Maintenance Manual

## *Document History*


<i>Document Number</i>	<i>BIO:POLICY:23</i>
<i>Created By</i>	<i>Kevin Corke</i>
<i>SUPERSEDED</i>	<i>Steven Haynes</i>
<i>Version</i>	<i>1.1</i>
<i>Date</i>	<i>19.02.2019</i>
<i>Review Date</i>	<i>19.02.2020</i>
<i>Approved by</i>	<i>Kevin Corke</i>




The  
University  
Of  
Sheffield.

## 1. Location

The facility is located on C road of the Royal Hallamshire hospital (RHH) site adjacent to the Biorepository freezer store.



### 1.1 Site Map



### 1.2 LN2 Store

## **2. Access & Security Systems**

The facility is maintained on the STH multi-lock system with a dedicated suited key. The number of keys is restricted and a signing in procedure is in place and only designated staff may enter the facility (see related SOP: Access).

The Facility is also monitored by CCTV which is monitored 24/7 by STH control room.

## **3. Heating and Ventilation**

The room has back ground heating supplied by Dimplex 1kW wall heaters to maintain an internal temperature above freezing.

Ventilation: the room is well ventilated with Louver doors and the addition of emergency extract fans which in the event of an Oxygen depletion alarm switch on automatically.

## **4. Fire prevention**

The Facility has smoke detectors linked to the Trust Building management system which is manned 24/7. Appropriate fire extinguishers are also located within the facility

## **5. Associated Documents**

	Document	Document Reference
1	LN2 freezer store O & M Jencons Manual	BIO:POLICY:24
2	LN2 Store security & Access SOP	BIO:SOP:30
3		
4		
5		
6		
7		